

# COMPOUND STUDENT MICROSCOPE BM-4

This microscope have well designed heavy and sturdy pressure die casted monocular body inclinable to 90° to meet the requirements of School and college students



**Eye Piece:** Wide-Field 10X Eye-piece or set of 10x & 15x Huygenian Eye-pieces.

**Nosepiece:** Triple revolving nosepiece with click stops.

**Stage:** High quality S.S. Stage clips provided to hold the slide.

**Focussing:** By vertical movement of the observation tube, with coarse adjustment using rack-pinion movement and fine motion by screw level mechanism.

**Condenser:** Abbe condenser having numerical aperture of 1.25 with iris diaphragm and swing out filter holder can be moved through rack and pinion. This condenser unit incorporates high efficiency optical system for optimum utilisation of light.

**Illumination:** Arc-Reflector set for using natural light to study the slide.

**Objective:** Achromatic 10x and 45x.

**Packaging:** Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.

