

COMPOUND STUDENT MICROSCOPE **BM-3**



This microscope have well designed heavy and sturdy pressure die casted monocular body inclinable to 90° to meet the requirements of School and college students

Eye Piece: Wide-Field 10X Eye-piece or set of 10x & 15x Huygenian Eye-pieces.

Nosepiece: Triple revolving nosepiece with click stops.

Stage: High quality S.S. Stage clips provided to hold the slide.

Focusing: By vertical movement of the observation tube, with coarse adjustment using rack-pinion movement and fine motion by screw level mechanism.

Condenser: Abbe condenser having Numerical aperture of 1.25 with iris diaphragm and swing out filter holder.

Illumination: Arc-Reflector set for using natural light to study the slide.

Objective: Achromatic 10x and 40x.

Packaging : In mould styrofoam box or wooden storing cabinet with lock and key (on demand)

Packaging: Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.

