



# MEDICAL MICROSCOPE BM-5



MONO



TL



ABBE



2 DAY

ISI Marked Medical Microscope has a well designed heavy and sturdy die casted monocular body inclinable to 90° to meet the needs of school and college students

- > Observation tube with mechanical tube length 160mm with triple nose piece.
- > Separate coarse and fine motion.
- > Fixed Co-axial horizontal graduated mechanical stage 120mm x 120mm with convenient adjustment for slide movement.
- > Double lens bright field abbe condenser NA 1.25 is fixed below the stage with iris diaphragm and swing out filter holder , can be moved through rack-pinion.
- > Adjustable down stopper for preventing accidental damages to slide
- > Illumination by piano concave mirror

**Also available without ISI mark**

**Packaging:** Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.

