

## Pathological Binocular Microscope

- > Two 10x Wide Field Eye -piece.
- > 4 nos. Objectives, Parafoal Achromatic, 4x, 10x, 40x SL and 100x SL Oil immersion.
- > Stage : Horizontal Mechanical of size 120 X 125mm min.
- > Illumination : LED Light with on/off switch, light intensity regulator and Plano concave mirror attachment suitable for 220V, 50 Hz. AC supply.
- > Eye Pieces for Binocular Microscope confirm to IS : 8275:1976

**Also available without ISI mark**

**Packaging:** Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.

# PATHOLOGICAL BINOCULAR MICROSCOPE BM-6bi

