

- > Ergonomic Design
- > Exceptional Precision
- > Accurate Measurements
- > Highest Industry Standard
- > Outstanding Optical Performance

The BM-9bi reflects friendly design suitable for both the classroom and laboratory environment. Embodied with sturdy C-shaped aluminum alloy die-cast stand, high quality optics and optical coatings, and brass gears, all of the which compose the BM-9bi to resist the toughest of user environments.

- > Binocular Head having Hard coated prisms
- > Quadruple revolving Ball Bearing Nose piece
- > Co-axial Focussing system
- > Double layered mechanical stage with low positioned coaxial knobs
- > Bright LED Light with Solid state power supply circuitry.
- > Semi-Plan Achromatic objectives, Extra Wide Field Eye Piece (Hard Coated) along with high quality condenser lens for better results for all magnifications.

Objectives: 4x, 10x, 40x, 100x (oil)

Eye pieces : EWF 10x (18mm Paired)

Also available without ISI mark

Packaging: Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.

BINOCULAR MICROSCOPE BM-9bi



BM-9TR