

Trinocular Gemological Microscope adopts streamline integrated mechanical design. Easy operation, long working distance, clear resolved image, enclosure integrated. It can be extensively used for observation studies of Gems and Stones.

Specifications:

- Superior Japanese (Hoyo or other) optics
- 20X-80X zoom stereo head
- True darkfield illumination halogen lamp 12V 10W
- Adjustable iris diaphragm
- Gem clip
- Eye piece : 10/20x
- Internationally adaptable voltage 120V/240V
- 360 degrees rotating head
- Dist cover, One year mechanical and electrical warranty

GEMOLOGICAL MICROSCOPE

GM-01

