



ALMICRO Inclined Metallurgical Microscope comes fitted with Monocular observation tube inclined at 45°.

It is provided with Bright field horizontal illuminator with centering provision, field and aperture diaphragm, and slots for dropping filter. It comes equipped with Quadruple revolving nose piece.

There is a Special graduated mechanical stage for easy study of the samples and is supplied with variable light control transformer and following optics in box :

LWD Objectives: M10X and M45X (SL)

Wide Field Eye pieces: WF10X and H15X

Packaging: Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.

Inclined Metallurgical Microscope

MM-2

