

Presenting our latest model of trinocular inverted metallurgical microscope having co-axial coarse and fine motion assembled on ball bearing guides. It is equipped with horizontal illuminator with centering provision, field and aperture diaphragm and slots for dropping filters. Built in base, variable solid state transformer. Large sturdy stage with X-Y axis movements for optimum study of sample.

Objectives: M5x, M10x, M20x, M40x/45x

Eye pieces: WF 10x (Paired)

Optional Objectives: 100x

Magnification : 50x to 450x

Packaging: Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.

Inverted Metallurgical Microscope

MM-5

