



Description:

MICROPROCESSOR pH METER (3 POINT CALIBRATION)

TECHNICAL SPECIFICATION:

pH	: 0 to 14.00 pH
mV	: -1500 to + 1500Mv
Display	: 3½ Digit Display
Modes	: pH-mV-temperature
Temp. indication	: available
Calibration	: 3 point with auto buffer
Resolution	: 0.01 pH
Repeatability	: + 0.01 pH ; + 1 mV
Accuracy	: pH : + 0.01 pH + digit
: mV	: + 1 mV + 1 digit
Temp. Compensation	: Automatic
Temp. Compensation Range	: 0 - 100 °C
Temp. Accuracy	: ±0.3
Operating humidity	: Up to 75%
Operating Temp.	: Up to 35°C
Display	: 3½ digit bright red seven segment LED display with automatic polarity and decimal point indication.
Power	: 230 V + 10 % AC , 50 Hz.
Dimensions	: 76 mm(H) X 275 mm (W) X 175 mm(D)

pH METER

Model : PH-10

