


- > Monocular head replaceable with Bertrand lens, Sliding Analyzer 1/2 & 1/4 wave compensator.
- > Co-axial Focussing system
- > Rotatable stage divided into 0-360° with vernier
- > Adjustable Condenser N.A. 1.25.
- > Polarizer duly graduated with in- out system. Built in light arrangement with 6V20W halogen lamp.
- > Solid State light control arrangement.
- > Paired eye piece 5x(cross), 10x(ocular), 15x (optional).
- > Objectives P5x, P10x, P40x/45x

Packaging: Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.



STUDENT POLARIZING MICROSCOPE PM-10