

INVERTED TISSUE CULTURE TM-8/8A

ALMICRO Inverted Tissue Culture Microscope is used for tissue culture examination. This Microscope has coarse and fine adjustment, incident LED Light with variable light control arrangement.

It comes equipped with long working distance objectives and binocular head, complete with following optics in box :

LWD Objectives: 5x, 10x, 20x, 40x

Wide Field Eye pieces: 10x (Paired)

Packaging: Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.



TM-8A